

CuroCell® A4 Cirrus Active Hybrid Air Mattress

Product Code: 36CCA4C

A hybrid mattress used as an aid in the prevention and treatment of pressure injuries up to and including stage 4 pressure injuries.



230kg
SWL



Very High Risk
Mattress

4

Aids in the Treatment of up
to & including **Stage 4** PI



Automatic
Pump



Anti Shear
Inner Cover

The CuroCell® A4 Cirrus Hybrid Mattress is used as an aid in the prevention and treatment of pressure injuries up to and including stage 4 pressure injuries. Specifically well suited for pain therapy, palliative or end of life individuals. The CuroCell® A4 Cirrus is constructed of air cells encapsulated by differing foams, including a soft foam surface that provides excellent user comfort and pressure redistribution. The integrated heel function/slope offers pressure redistribution on the heels, a known anatomical location prone to pressure injuries. The CuroCell® A4 Cirrus air mattress has been developed to allow healthcare professionals to focus on patient-centered care. Using AI and high-precision sensors, the system continuously monitors the individual's weight, height, and position, adapting to each individual automatically with no manual adjustment.

The Pulsating Mode™ uses artificial intelligence to combine soft, alternating movements with constant low pressure to offer a large contact area and reduce high peak pressures. Its clinical effectiveness is well-documented. The Pulsating Mode™ is developed to improve comfort and stability by reducing movements and fluctuations in the support surface. In addition, the silent running of the control units maintains a restful environment for undisturbed sleep and recuperation.

Product Features:

- Integrated heel function, reducing pressure on the heels
- Shear reduction and hygiene cover, the surface has a removable and liquid-proof hygiene cover manufactured with a 4-way stretch material to reduce the risk of shear forces on the skin. Additionally, the cover is vapour permeable and has waterfall flaps designed with infection control in mind
- CPR and transport function
- Smart cable management & easy handling
- Inner anti-shear cover
- Safe Working Load 230kg
- Quality product from Europe

Pump Features:

- Automatic AI-driven weight and position detection requiring no manual weight entry
- Three program modes: Pulsating, Alternating (GALP) and CLP (constant low pressure)
- Maintains air during transportation and power failure - the support surface air pressure can be maintained for at least 12 hours without being connected to the control unit
- Silent and vibration-free. The system reuses air, enabling the control unit to operate quietly and intermittently, supporting restful sleep and recovery
- Automatic 20-minute timer after max inflate/care mode is enabled
- WARRANTY: 2 Years



Specifications:

Product Code	3674A4	3675A4	3687A4	3676A4	3677A4	3846-1	36770A4
Size	Single	Single	King Single	King Single	Extra Wide	Double	Queen
Width	850 mm	900 mm	1000 mm	1050 mm	1200 mm	1350 mm	1500 mm
Length	2000 mm	2000 mm	2000 mm	2000 mm	2000 mm	2200 mm	2000 mm
Thickness	210 mm	210 mm	210 mm	210 mm	210 mm	210 mm	210 mm
SWL	230 kg	230 kg	230 kg	230 kg	230 kg	230 kg	230 kg

Issue Date: 240325
Revision Date: 290726 | V3