

CuroCell A4 CX16 Active Air Mattress

Product Code: 36CCCX16

The A4 CX16 Active Air Mattress is a fully automatic & flexible alternating pressure mattress replacement system focusing on user-centered care.



The CuroCell® A4 CX16 provides simplicity together with user comfort while effectively aiding in the prevention and treatment of pressure injuries up to and including category 4, and unclassified pressure injuries. It is designed for simplicity and provides features for excellent user comfort, reduction of shear forces and improvement of the microclimate.

The CuroCell® A4 CX16 mattress system includes an intelligent pump with a user-friendly interface. It is fully automatic without need for manual settings. The system identifies the user and automatically manages all settings to optimise the pressure according to the user's needs.



200kg
SWL



High Risk
Mattress

4

Suitable for
Category 4

AUTO

Product Features:

- Individually replaceable cells with safety mattress in HR-foam
- Heel function
- Multifunctional inner cover
- Unique CPR-function with transport function
- Smart cable management
- Unique design for easy handling and usage
- Well protected zipper
- Fully automated pump that adapts individually
- User friendly interface
- Audible and visual notifications
- Three program modes: Pulsating, Alternating and CLP (constant low pressure)
- Caring mode (max-firm)
- Pack&Go™ function
- Silent running and vibration free
- Low voltage (12V power supply)
- Soft friction hangers with extra soft tip
- Air flow control; Constant temperature, increased user comfort, energy reduction



Specifications:

Product Code	3626	3627	3628
Size	Single	Wide	King Single
Width	850 mm	1000 mm	1050 mm
Length	2000 mm	2000 mm	2000 mm
Thickness	160 mm	160 mm	160 mm
SWL	200 kg	200 kg	200 kg