

# CuroCell SAM PRO Static Air Mattress

Product Code: 36CCSMP

The CuroCell SAM PRO Static Air Mattress is a mattress Replacement System with very high comfort.



The CuroCell SAM PRO Static Air Mattress is a mattress replacement system with an integrated foam core can be used as an effective aid in the prevention and in the treatment of pressure injuries up to and including category 3. It can also be used as an aid for pain therapy.

This constant low pressure mattress offers an exceptionally high level of comfort as well as excellent pressure redistribution. With the vastly improved higher maximum weight capacity of 250 kg, this overlay will meet the ever increasing challenges of heavier patients.

The CuroCell S.A.M.® PRO features a new cost effective pump that includes both CPR and Transport functions. The pump automatically determines when to change pressure in contrast to running constantly.



**250kg**  
SWL



**High Risk**  
Mattress

**3**

Suitable for  
**Category 3**

**P**

**Prevention**  
Mattress

### Product Features:

- CPR function – rapid deflation
- Soft, elastic, vapour permeable hygienic cover
- Very good cleaning properties
- Silent running and durable pump
- Easy to handle and maintain
- Energy efficient
- Maximum Pressure (caring mode)
- Extra soft inner cover
- Heel function
- Power Supply of 230V / 50Hz
- Prevention and treatment up to and including category 3
- Warranty: 2 years
- Cleaning of cover: Wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying



### Specifications:

Product Code	3683	3684	3685
Width	850 mm	1050 mm	1000 mm
Length	2000 mm	2000 mm	2000 mm
Thickness	160 mm	160 mm	160 mm
SWL	250 kg	250 kg	250 kg