

CuroCell® SAM PRO Reactive Air Mattress

Product Code: 36CCSBP

The CuroCell® SAM PRO Reactive Air Mattress is a mattress Replacement System with very high comfort.



250kg
SWL



High Risk
Mattress

3

Aids in the Treatment of up
to & including **Stage 3** PI



Silent & Vibration
Free Pump



4-Way
Stretch Cover

The CuroCell® SAM PRO Reactive Air Mattress is a static air replacement mattress which has been proven to provide exceptional comfort, pressure and pain relief. It is used preventively and in the treatment of pressure injuries up to and including category 3. The CuroCell® SAM PRO combines the comfort of the foam with the pressure-relieving immersion and envelopment properties from the central air chamber.

The mattress is divided into three zones, which, amongst other advantages, give better support when in the seated position. The pump, together with the mattress, is easy to use and very quiet. The 20-minute automatic timer embedded in the pumps max inflate function ensures safety surrounding immersion and envelopment principles when utilising this function for cares and transfers. Designed with an anti-shear inner cover and a soft, vapour-permeable outer cover, the CuroCell® SAM PRO mattress has very well-documented hygiene standards.

Product Features:

- Suitable for high-risk individuals
- Reactive (static) air mattress replacement system
- Air cells are produced from soft, durable polyurethane, not PVC
- 20-minute integrated timer for max inflate mode
- Vapour permeable, water resistant and 4-way stretch outer cover
- Significant noise reduction when cells move compared to other market products
- Mattress comes complete with an anti-shear inner cover, which has been proven to reduce pressure by up to 8%
- Easy-to-use dial pump with a table to calculate the most appropriate setting
- Safe Working Load 250kg
- Quality product from Europe
- WARRANTY: 2 years



Specifications:

Product Code	3683	3684	3685
Width	850 mm	1050 mm	1000 mm
Length	2000 mm	2000 mm	2000 mm
Thickness	160 mm	160 mm	160 mm
SWL	250 kg	250 kg	250 kg

Issue Date: 171224
Revision Date: 290726 | V3