

Staebel Zonecare 80 Medical Mattress

Product Code: 372SZC

The Staebel Zonecare 80 mattress is a high quality foam mattress featuring a profile-cut pressure relieving foam insert.



150kg
SWL



Low Risk
Mattress



Prevention
Mattress



Firm
Surrounds



Multi-Way
Stretch Cover

The Staebel Zonecare 80 Medical Mattress is a quality foam mattress featuring a profile-cut/castellated pressure-relieving foam surface, enhancing immersion and envelopment for the user, and providing pressure redistribution. The firm surround design additionally supports independent bed mobility, including transfers. The Staebel Zonecare 80 Medical Mattress is considered a comfortable yet cost-efficient solution for both home care and residential aged care settings and is a reversible support surface suitable for low-pressure injury risk individuals.

Mattress Features:

- Castellated foam surface to enhance envelopment and provide pressure redistribution
- Suitable for low-risk individuals
- Firm surrounds to enable independent bed transfers and mobility

Cover Features:

- Waterfall zipper cover to minimise fluid penetration and infection control risks
- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Water-resistant and fire-retardant cover
- Multi-way stretch cover
- Mattress cover has welded seams to meet strict infection control standards
- WARRANTY: Foam core – 5 years , Cover – 1 year



Specifications:

Product Code	3728Z	3729Z-NINO	3729Z-W	3728Z-KS	372SZC-K
Size	Single	Nino	Wide	King Single	King Double
Width	850 mm	900 mm	950 mm	1050 mm	1800 mm
Length	1980 mm	1970 mm	1980 mm	1980 mm	2000 mm
Thickness	150 mm	150 mm	150 mm	150 mm	150 mm
SWL	150 kg	150 kg	150 kg	150 kg	150 kg

Issue Date: 120224
Revision Date: 020326 | V3