

# Optimal 5zon Pressure Care Mattress with Evac Handles

Product Code: 37E5Z

A flexible mattress replacement system for pressure injury prevention using our patented DuoCore® and a two-part foam core.



The Optimal 5zon with Evac Handles is intended for all healthcare environments as an aid in the prevention and treatment of pressure injuries up to and including category 2. It is a flexible mattress replacement that is placed directly on the bed and replaces the standard mattress and mattress pad.

Optimal 5zon® has good pressure redistribution properties and is also equipped with our patented DuoCore® sit function. The sit function has been designed to provide support when seated, while also offering good comfort in lying position.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable, which contributes to reduce the risk of skin maceration.



180kg  
SWL



High Risk  
Mattress

2


Suitable for  
Category 2

P

Prevention  
Mattress

### Product Features:

- Elastic material and a zone and module system gives pressure equalisation.
- Function: as an aid for prevention and treatment up to and including category 2
- The mattress is specially designed to reduce the effects of shear forces.
- Our patented sit function DuoCore® is designed to give support when sitting and also to provide good comfort when lying.
- The hygienic cover is vapour permeable, which contributes to reduce the risk of skin maceration.
- The cover is easy to clean and has a good hygienic standards.
- Machine wash max 95 °C, tumble drying
- Safe Working Load 180kg



### Specifications:

Product Code	3770	3772	3774	3776
Width	850 mm	900 mm	1050 mm	1200 mm
Length	2000 mm	2000 mm	2000 mm	2000 mm
Thickness	140 mm	140 mm	140 mm	140 mm
SWL	180 kg	180 kg	180 kg	180 kg