

Staebel Gel Segment Bed Pad

Product Code: 3853

The Staebel Gel bed segments are made of Akton polymer which mimics human fatty tissue in consistency; providing equally distributed pressure management.



Discover superior comfort and pressure relief with Stabeel Gel Segment Bed Pads, designed to enhance your sleeping experience and support personalised care needs. Made from innovative Akton polymer, these overlays replicate the consistency of human fatty tissue, delivering consistent and effective pressure management without the migration issues commonly found in liquid gel alternatives.

These overlays provide flexibility to target specific areas like shoulders, hips, or feet, or to cover the full surface of a standard mattress. Use one, two, or all three segments to suit your unique requirements, ensuring maximum support where it's needed.



**Unlimited
SWL**



**Pressure
Distribution**



**Bariatric
Needs**



**Prevention
Product**

Product Features:

- **Premium Akton Polymer:** Mimics human fatty tissue, offers uniform pressure distribution, and resists migration for reliable support.
- **Customisable Coverage:** Use individually or in combination for targeted pressure relief or full-body coverage.
- **Flexible Sizing:** Dimensions of each segment are 69 x 63 x 1.8 cm, with three segments covering the full length of a single bed.
- **Easy Maintenance:** Durable and simple to clean, requiring no regular upkeep.
- **Weight Compatibility:** Supports users of all weights without restriction, ensuring universal usability.
- **Lightweight & Manageable:** Each segment weighs 8kg, allowing for easy handling and repositioning.

Specifications:

Product Code	3853
Width	630 mm
Length	690 mm
Thickness	18 mm
SWL	Unlimited

