

Staebl Zonecare Mental Health Mattress

Product Code: 3M9000MH

Specifically engineered for use in high-security mental health environments, where durability, safety, and hygiene are essential.



Mattress Only



200kg
SWL



Mental Health
Use



Prevention
Mattress



Highest Infection
Control

Constructed with a single-core foam interior, the mattress supports four-way turning to extend product life and promote even wear. Its heavy-duty PVC cover is fully sealed and enclosed, providing a robust barrier that protects the foam core from damage and contamination.

The cover features a double-laminated surface to regulate internal pressure regulation and ensures mattress longevity. High-frequency welded seams ensure the mattress is waterproof, preventing fluid ingress and supporting strict infection control protocols. Designed to withstand challenging conditions, the Zonecare Mental Health Mattress is hypoallergenic, anti-static, fire retardant, and resistant to tearing, stains, and microbial growth, making it a reliable choice for secure healthcare settings.

Product Features:

- High-Security Construction: Fully sealed, enclosed cover with high-frequency welded seams to prevent tampering and protect the foam core from willful damage
- Enhanced Infection Control: 100% fluid-proof, anti-microbial, stain-resistant surface that meets and exceeds hygiene requirements
- Durability and Longevity: Single-core foam with four-way turn capability and integrated air valve for extended use and structural integrity
- Safety-Compliant Materials: Fire-retardant, hypoallergenic, and anti-static properties for patient and staff safety
- Safe Working Load 200kg
- WARRANTY: Foam core – 2 years , Cover – 2 years



Specifications:

Produce Code	3M9000MH-S	3M9000MH-KS	3M9000MH-EW
Size	Single	King Single	Extra Wide
Overall Width	860 mm	1050 mm	1200 mm
Overall Length	1980 mm	1980 mm	1980 mm
Thickness	150 mm	150 mm	150 mm
SWL	200 kg	200 kg	200 kg

Issue Date: 231224
Revision Date: 190326 | V2