

Staebel Zonecare 40 Medical Mattress

Product Code: 3M9040

The Staebel Zonecare 40 Medical Mattress features a profiled foam surface designed to enhance the free flow of air and reduce surface tension. It comes with a two-way stretch vapour permeable cover that has sewn seams and a zipped closure, ensuring durability and ease of maintenance.



The Staebel Zonecare 40 Medical Mattress features a profiled foam surface designed to enhance the free flow of air and reduce surface tension. It comes with a two-way stretch vapour permeable cover that has sewn seams and a zipped closure, ensuring durability and ease of maintenance. The mattress is profiled on both sides for added comfort and support.

Product Features:

- Profiled foam surface
- 2 Way stretch cover
- Vapour permeable cover
- Sewn Seams
- 4 way turn
- Zipped closure, profiled on both sides
- Flame resistant cover and foam
- Range of sizes available



Specifications:

Product Code	3M9040-EN	3M9040-N	3M9040-S	3M9040-KS	3M9040-EW	3M9040-D	3M9040-Q	3M9040-K
Size	Extra Narrow	Narrow	Single	King Single	Extra Wide	Double	Queen Double	King Double
Width	670 mm	720 mm	870 mm	1070 mm	1150 mm	1370 mm	1520 mm	1780 mm
Length	1960 mm	1960 mm	1960 mm	1960 mm	1960 mm	1960 mm	1960 mm	1960 mm
Thickness	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm
SWL	150 kg	150 kg	40 - 130 kg	40 - 130 kg	150 kg	40 - 130 kg	40 - 130 kg	150 kg