

Staebel Zonecare 40 Medical Mattress

Product Code: 3M9040

The Staebel Zonecare 40 Medical Mattress features a profiled surface that promotes the free flow of air and reduces surface tension, designed to deliver a cost-effective solution across healthcare environments.



140kg
SWL



Low Risk
Mattress



Prevention
Mattress



Profiled
Foam Surface



Multi-Way
Stretch Cover

The Staebel Zonecare 40 Medical Mattress offers a two-way stretch cover with sewn seams and a zipped closure. Additionally, the cover is vapour-permeable to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management. The mattress is rotatable, available in a range of sizes and has been designed to deliver a cost-effective solution across healthcare environments.

Mattress Features:

- Suitable for low-risk individuals
- Profiled foam surface to support envelopment and therefore pressure redistribution
- 4-way turn
- Safe Working Load of 140kg

Cover Features:

- Multi-way stretch cover
- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Waterfall zipper cover to minimise fluid penetration and infection control risks
- Flame-resistant cover
- WARRANTY: Foam core - 5 years , Cover - 1 years



Specifications:

Product Code	3M9040-EN	3M9040-N	3M9040-S	3M9040-KS	3M9040-EW	3M9040-D	3M9040-Q	3M9040-K
Size	Extra Narrow	Narrow	Single	King Single	Extra Wide	Double	Queen Double	King Double
Width	670 mm	720 mm	870 mm	1050 mm	1150 mm	1350 mm	1480 mm	1780 mm
Length	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm
Thickness	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm
SWL	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg

Issue Date: 231224
Revision Date: 190326 | V2