

Staebel Zonecare 50 Medical Mattress

Product Code: 3M9050

The Staebel Zonecare 50 Medical Mattress consists of three layers of pressure-sensitive memory foam to provide superior pressure care designed to deliver a cost-effective solution across healthcare environments.



50-150kg
SWL



Low Risk
Mattress



Prevention
Mattress



Firm
Surrounds



Multi-Way
Stretch Cover

The Staebel Zonecare 50 Medical Mattress consists of three layers of pressure-sensitive memory foam at differing densities to provide pressure re-distribution. The lower layer is a high-density foam construct, so the patient won't "bottom out." The mattress can accommodate a wide range of user weights, and the firm surrounds enable independent bed transfers and mobility. The cover is durable, water-resistant and vapour-permeable to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management. Lastly, the base of the is made of durable PVC to promote longevity and prevent product wear due to sliding. This mattress is available in a range of sizes and has been designed to deliver a cost-effective solution across healthcare environments.

Mattress Features:

- Suitable for low-risk individuals
- Firm surrounds to enable independent bed transfers and mobility
- 3-layer graded foam structure
- Visco pressure-sensitive memory foam top layer
- Resilient mid-core
- High-density foam lower layer to prevent bottoming out
- Effective pressure redistribution
- SWL 50-150kg

Cover Features:

- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Heavy-duty PVC base
- Multi-way stretch cover
- Waterfall zipper cover to minimise fluid penetration and infection control risks
- WARRANTY: Foam core - 5 years , Cover - 2 years



Specifications:

Product Code	3M9050-EN	3M9050-N	3M9050-NINO	3M9050-S	3M9050-KS	3M9050-EW	3M9050-D	3M9050-Q	3M9050-K
Size	Extra Narrow	Narrow	Nino Cot	Single	King Single	Extra Wide	Double	Queen	King Double
Width	670 mm	720 mm	900 mm	860 mm	1050 mm	1150 mm	1350 mm	1480 mm	1780 mm
Length	1980 mm	1980 mm	1970 mm	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm	1980 mm
Thickness	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm	150 mm
SWL	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg	50 - 150 kg