

# Staebel Zonecare 50 Medical Mattress

Product Code: 3M9050

The Staebel Zonecare 50 Medical Mattress consists of three layers of pressure-sensitive memory foam to provide superior pressure care designed to deliver a cost-effective solution across healthcare environments.



50-150kg  
SWL



Low Risk  
Mattress



Prevention  
Mattress



Firm  
Surrounds



Multi-Way  
Stretch Cover

The Staebel Zonecare 50 Medical Mattress consists of three layers of pressure-sensitive memory foam at differing densities to provide pressure re-distribution. The lower layer is a high-density foam construct, so the patient won't "bottom out." The mattress can accommodate a wide range of user weights, and the firm surrounds enable independent bed transfers and mobility. The cover is durable, water-resistant and vapour-permeable to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management. Lastly, the base of the cover is made of durable PVC to promote longevity and prevent product wear due to sliding. This mattress is available in a range of sizes and has been designed to deliver a cost-effective solution across healthcare environments.

Mattress Features:

- Suitable for low-risk individuals
- Firm surrounds to enable independent bed transfers and mobility
- 3-layer graded foam structure
- Visco pressure-sensitive memory foam top layer
- Resilient mid-core
- High-density foam lower layer to prevent bottoming out
- Effective pressure redistribution
- SWL 50-150kg

Cover Features:

- Vapour Permeable cover to minimise the negative impact of moisture accumulation at the skin and surface interface, supporting microclimate management
- Heavy-duty PVC base
- Multi-way stretch cover
- Waterfall zipper cover to minimise fluid penetration and infection control risks
- WARRANTY: Foam core - 5 years , Cover - 1 years



Specifications:

| Product Code | 3M9050-EN    | 3M9050-N    | 3M9050-NINO | 3M9050-S    | 3M9050-KS   | 3M9050-EW   | 3M9050-D    | 3M9050-Q    | 3M9050-K    |
|--------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Size         | Extra Narrow | Narrow      | Nino Cot    | Single      | King Single | Extra Wide  | Double      | Queen       | King Double |
| Width        | 670 mm       | 720 mm      | 900 mm      | 870 mm      | 1050 mm     | 1150 mm     | 1350 mm     | 1480 mm     | 1780 mm     |
| Length       | 1980 mm      | 1980 mm     | 1970 mm     | 1980 mm     | 1980 mm     | 1980 mm     | 1980 mm     | 1980 mm     | 1980 mm     |
| Thickness    | 150 mm       | 150 mm      | 150 mm      | 150 mm      | 150 mm      | 150 mm      | 150 mm      | 150 mm      | 150 mm      |
| SWL          | 50 - 150 kg  | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg | 50 - 150 kg |

Issue Date: 231224  
Revision Date: 190326 | V2