

Staebel Zonecare 500 Bariatric Medical Mattress

Product Code: 3M9500


The Staebel Zonecare 500 Bariatric Medical Mattress provides care, pressure relief and comfort for bariatric users weighing up to 500kg. This mattress is the fitting choice for clinicians and facilities looking to achieve comfort, longevity and effective user care.



The Staebel Zonecare 500 Bariatric Medical Mattress provides care, pressure relief and comfort for bariatric patients weighing up to 500kg. Through the use of proven foam combinations, unique construction methods and high performance cover materials, this mattress is the fitting choice for clinicians and facilities looking to achieve comfort, longevity and effective patient care.

Product Features:

- 4-Layered Graded Pressure Care Core
- High specification foam for longevity and performance
- Provides graded postural support and patient comfort
- Foam conforms to AS/NZS 4088.1 (Int):1993 and FMVSS302 & SAE T2002
- Anti-microbial treatment applied to foam maximises longevity and core micro-climate



Specifications:

Size	King Single	Extra Wide	Double	Queen Double
Width	1050 mm	1150 mm	1350 mm	1480 mm
Length	1960 mm	1960 mm	1960 mm	1960 mm
Thickness	175 mm	175 mm	175 mm	175 mm
SWL	500 kg	500 kg	500 kg	500 kg