

Staebel Zonecare 500 Bariatric Medical Mattress

Product Code: 3M9500

The Staebel Zonecare 500 Bariatric Medical Mattress is specifically designed to support bariatric patients weighing up to 500kg, offering pressure relief, postural support, and long-term durability.



500kg
SWL



Low Risk
Mattress



Prevention
Mattress



Bariatric
Use



4 Layered
Foam Construct


Built with a four-layered graded foam core, the Staebel Zonecare 500 Bariatric Medical Mattress combines a high weight capacity with patient comfort, making it ideal for high-dependency and long-term care settings. Each layer works in unison to promote pressure redistribution, reduce shear, and enhance stability.

Using proven foam formulations, the Zonecare 500 delivers consistent performance under high load conditions. The foam core conforms to key safety and quality standards (AS/NZS 4088.1 (Int):1993 and FMVSS302 & SAE T2002) and includes anti-microbial treatment to protect both patient and product longevity. Paired with high-performance cover materials, this mattress meets the rigorous demands of clinical use while prioritising patient well-being.



Product Features:

- 4-Layered graded foam core
- Provides postural support and patient comfort
- Foam conforms to AS/NZS 4088.1 (Int):1993 and FMVSS302 & SAE T2002
- Anti-microbial treatment applied to foam maximises longevity
- Safe Working Load 500kg
- WARRANTY: Foam core - 5 years , Cover - 2 years



Specifications:

Product Code	3M9500-105	3M9500-115	3M9500-135	3M9500-148
Size	King Single	Extra Wide	Double	Queen Double
Width	1050 mm	1150 mm	1350 mm	1480 mm
Length	1980 mm	1980 mm	1980 mm	1980 mm
Thickness	200 mm	200 mm	200 mm	200 mm
SWL	500 kg	500 kg	500 kg	500 kg

Issue Date: 231224
Revision Date: 020326 | V3