

Sertain Hi-Lo Care Chair

Product Code: 5142

Versatile in design and function, this chair is a popular choice for many wards including intensive care units.



Versatile in design and function, our Sertain HILO Care Chair is a popular choice for many wards including intensive care units. The HILO function of the chair is a nurse's dream, allowing you to adjust the height of the chair to a position comfortable to assist the patient.

Product Features:

- Electric hi-lo ; low enough to stand from with ease, but high enough for a comfortable working height for carers
- Electric reclining backrest – full recline for supine position
- Electric tilting seat and elevating legrest
- Handset operated with preset auto chair and auto layflat functions
- Pressure management seating – gel cushion in seat, air in backrest and in seat, air in backrest
- Fold out, extendable footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Removable two-position armrests
- Swingaway three position headwings
- Adjustable head rest cushion
- Battery back up system
- Quality Australian made product
- Optional Colours available (see below)



Specifications:

Product Code	5142	5145	5147
Size	Standard	Large	Bariatric
Seat Width	520 mm	620 mm	720 mm
Seat Depth	520 mm	570 mm	600 mm
Seat Height	540 mm	550 mm	550 mm
Overall Height	550 - 850 mm	550 - 850 mm	550 - 850 mm
SWL	160 kg	300 kg	350 kg