

Regency Care Chair

Product Code: 5050G

The ultimate care chair designed to meet the needs of both the carer & the individual.



The Regency Care Chair is designed to give the user comfort by combining the original tilt system with an independent multi-positional back rest. Its many features take into account the needs of both carer and the individual of the chair. The Regency derives its status from a proven and successful track record across a broad spectrum of the industry and has features such as:

- Forward transfer friendly
- Maneuverable castor base
- High quality upholstery
- Side Transfer Friendly
- Superior range of tilt and completely flat recliner

Helping You Choose the Right Care Chair

Not all care chairs are designed for the same purpose. This guide supports clinical and procurement teams in selecting seating that aligns with the individuals needs, staff workflows, and long-term care planning. Points to consider when choosing a care chair include:



Level of postural and functional support required



Transfer methods (standing, lateral, hoist-assisted)



Clinical workload and monitoring requirements



Infection control and ease of cleaning



Long-term durability and maintenance



Resident dignity and environment aesthetics



When is the Regency Care Chair the Right Choice?

Regency Care Chairs are ideal when the individual requires:

- High or changing postural support
- Upright functional seating
- Standing or lateral transfer capability
- Reduced falls risk
- Integration into physiotherapy or mobility programs

Clinical Note: Pressure Monitoring & Air Seating:

Air-based seating systems may require regular clinical checks to ensure correct pressure or immersion depth. Without consistent monitoring, risks may include:

- Bottoming out
- Increased shear and friction
- Loss of positioning integrity
- Additional staff workload

Regency's structured foam and contoured support system provides consistent seating performance without ongoing adjustment.

Product Features:

- Reclining backrest – full recline to supine position
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests
- Forward seat tilt for ease to stand, legrest angles under seat
- Elevating contoured legrest complete with fold out footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Safe Working Load 180kg
- Patient lifter friendly castor base
- 5 year limited warranty
- 175mm castors, braking on rear, directional lock on front
- Quality Australian made construction

Accessories – Added after time of order

- Standard meal/activity tray with PVC edge (Only available with Standard Armrest)
- Oxygen Bottle and IV Kit
- CPR Layflat Bar

Return on Investment Over Time

While some chairs may appear cost-effective initially, long-term value should account for:

- Product lifespan
- Cushion replacement frequency
- Maintenance requirements
- Staff training and compliance effort



Swing-Away Armrests / Wings
enhances side transfers



Layflat Position
Function



2-Position Armrests
for upright posture and support



Suitable for
Pressure Management

Comparison Chart:

User Requirements	Regency Care Chair	Air Chair
Bariatric Support	✓	Model dependent
High Falls Risk	✓	Limited
Standing Transfers	✓	✗
Palliative Care	Case dependent	✓
Structured Posture	✓	✗
Reduced Monitoring Workload	✓	✗
Primary Purpose	Structured posture & functional seating	Comfort & pressure care
Seating System	Contoured foam, gel, air backrest	Air or immersion layers
Postural Support	Excellent	Limited
Tilt & Recline	Full tilt & flat recline	Recline (model dependent)
Leg & Footplate Access	Flip-up for transfers	Often fixed
Drop-Down Arms	✓	Usually no
Standing Transfer Friendly	✓	✗
Pressure Monitoring	Not required	Required
Durability	High	Moderate
Infection Control	Excellent	Moderate
Aesthetics	Furniture Like	Institutional

We are
activeMEDICAL

P. 1800 267 267

E. hello@activemedical.au

www.activemedical.au