

Regency R2900 Care Chair

Product Code: 5040

This Regency R2900 Care Chair is designed to be carer and client friendly.





This Regency R2900 Care Chair is designed to be carer and client friendly by incorporating the original tilt seat system into a multi-positioning, multi-functional and manoeuvrable care chair. The Regency R2900 derives its status from a proven and successful track record across a broad spectrum of the industry.

Product Benefits:

- Reclining backrest – full recline to supine position
- Tilting seat – for ease to stand, legrest angles under seat
- Pressure management seating – gel cushion in seat, self-air adjusting back, profiled head cushion
- 2-position armrests
- Elevating contoured legrest complete with fold out, extendible footplate
- Low skin shear 2-way stretch fabric on seating surfaces
- Drop down armrests and swingaway headwings – side transfer friendly
- Covers replaceable
- Patient lifter friendly castor base
- Australian Made & Owned
- Safe Working Load of 180kg

Specifications:

Product Code	5040
Seat Width	460 mm
Seat Depth	540 mm
Seat Height (Tilt/Normal)	520 mm / 550 mm
Backrest Height	680 mm
Armrest Height	175 mm
Legrest Length	470 – 550 mm
Overall Width	660 mm
Overall Depth	1050 mm
Overall Height	1280 mm
SWL	180 kg