


Shoalhaven AirCell Recliner

Product Code: 5230

Ensuring maximum comfort for those confined to a chair, whilst allowing the carer to operate the chair with maximum safety and care.





Following on the success of the Shoalhaven comes the Shoalhaven Aircell. This has been engineered to ensure maximum comfort for those confined to a chair, whilst allowing the carer to operate the chair with maximum safety and care. Soft, stretch fabric, quality pressure sensitive foam and a shaped backrest giving lumbar and lateral support provide ultimate comfort.

Product Benefits:

- Air cell cushioning combined with slow release memory foam, for excellent pressure care
- Compact postural design with tilt in space recline action, and high legrest elevation for reduced fluid retention
- Multiple stretch, water proof fabric, offering low skin shear, anti bacterial and non allergenic properties
- Designed for use with standing lifters & transfer hoists
- Hospital quality 100mm locking castors
- Head rest cushion included
- Manufactured in Australia

Specifications:

Product Code	5230
Seat Width	520 mm
Seat Depth	500 mm
Seat Height	500 mm
Backrest Height	750 mm
Arm Height	220 mm
Legrest Length	400 mm
Overall Width	780 mm
Overall Depth	1000 mm
Overall Height	1170 mm
SWL	130 kg