

# Verity Mobile Recliner

Product Code: 5266

Offering healthcare-grade durability and exudes an inviting, homelike aesthetic that brings warmth to clinical environments.



The Verity Recliner offers an inviting, homelike aesthetic that brings a touch of residential appeal to clinical environments. Designed with both comfort and functionality in mind, it features a three-position recline with an easy-to-operate mechanism, allowing users to adjust their seating with minimal effort. Ergonomically-positioned push handles make transport simple and efficient, enhancing mobility for caregivers. Additionally, the recliner’s stitching method is specially designed to reduce seams, helping to increase infection resistance and support a cleaner, more hygienic setting.

**Product Benefits:**

- 4 locking castors
- 3 position recline with easy to operate mechanism
- Ergonomically-positioned push handles for patient transport
- Stitching method reduces seams for increased infection resistance
- CAL TB133 compliance
- Moisture Barrier
- Optional arm caps available in polyurethane
- Contrasting fabrics are obtainable

**Optional Features:**

- Optional wood, polyurethane, or solid surface arm caps
- Choice of locking casters, wood legs or metal legs
- Side-mount tablet with 360° swivel
- Flip-up side tablet
- Custom finishes
- Contrasting fabrics
- Moisture barrier

**Specifications:**

|             |         |
|-------------|---------|
| Width       | 740 mm  |
| Depth       | 740 mm  |
| Height      | 1120 mm |
| Seat Height | 480 mm  |
| Seat Width  | 585 mm  |
| Arm Height  | 660 mm  |
| SWL         | 225 kg  |

