

Bowen Lever Recliner

Product Code: 5513

The Bowen range of chairs have been adopted for hospital use due to their comfortable seating and range of options.



The Bowen range of chairs have been adopted for hospital use due to their comfortable seating and range of options. Economically priced, the Bowen features a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise.

Product Benefits:

- Single piece vinyl backrest for infection control
- Upholstered legrest with calf support
- No sag springs in seat and lumbar for long lasting support
- Tasmanian Oak hardwood frame for durability
- High quality Dunlop Enduro foams used
- Trigger Release Function
- Quality Australian made product
- Upholstered in Austex Studio Encore Vinyl range
- Available with Mobile Base

Specifications:

Product Code	5513
Seat Width	480 mm
Seat Depth	530 mm
Seat Height	480 mm
Backrest Height	600 mm
Overall Width	880 mm
Overall Depth	810 mm
Overall Height	1070 mm
SWL	120 kg