

Ashley Luxor Optima Manual Recliner

Product Code: 5519ALO

The Ashley Luxor range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated.



The Ashley Luxor range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated. This is achieved with layered pressure reducing foam and stretch PU upholstered on the seating surfaces. All the favourite features of the standard Ashley are still included – a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. Australian made and available in manual lever recliners and also electric lift/recline chair models in one or two motor models.

Product Benefits:

- Upholstered leg rest with calf support
- Manual operated recliner action with lever
- No sag springs in seat and lumbar for long lasting support
- Quality Australian made
- Waterfall backrest with Zips
- 4x Lockable castors
- Multi layer foam pressure care seating
- Combination of 2-way stretch upholstery
- Warranty: Fabric – 12 months, Frame – 24 months



Specifications:

Product Code	5519MI	5519	5519MX
Size	Mini	Medium	Maxi
Seat Width	480 mm	480 mm	530 mm
Seat Depth	490 mm	500 mm	520 mm
Seat Height	540 mm	560 mm	580 mm
Backrest Height	600 mm	640 mm	680 mm
Overall Width	780 mm	780 mm	830 mm
Overall Height	1100 mm	1130 mm	1230 mm
SWL	120 kg	120 kg	120 kg

Allow +/- 20mm tolerances