

Ashley Luxor Optima Recliner, One Motor

Product Code: 5522ALO

The Ashley Luxor Optima range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated.



The Ashley Luxor Optima range of chairs have been specially developed for those with pressure care needs or those needing increased comfort when seated. This is achieved with layered pressure reducing foam and stretch PU upholstered on the seating surfaces. All the standard Ashley are still included – a waterfall style backrest, sculptured armrests, side pockets and padded leg chaise. The two motor action enables independent control of backrest and leg rest operation for ultimate flexibility of chair position.

Product Benefits:

- Upholstered leg rest with calf support
- One motor wall saver action – chair can be positioned close to within 100mm from wall and still recline
- Backrest reclines to half lying down position
- Chair raises to assist standing transfer
- No sag springs in seat and lumbar for long lasting support
- Tasmanian Oak hardwood frame for durability
- High quality foam used
- Battery back up capability
- Small castors enable chair to be wheeled (not patient transport)
- Quality Australian made product
- One Motor Design
- 2 year electronic warranty, 5 year frame warranty



Specifications:

Product Code	5522P	5522MI	5522	5522MX
Size	Petite	Mini	Medium	Maxi
Seat Width	440 mm	480 mm	480 mm	530 mm
Seat Depth	440 mm	490 mm	500 mm	520 mm
Seat Height	470 mm	470 mm	500 mm	530 mm
Backrest Height	570 mm	600 mm	640 mm	680 mm
Overall Width	780 mm	780 mm	780 mm	830 mm
Overall Height	1020 mm	1020 mm	1070 mm	1170 mm
SWL	120 kg	120 kg	120 kg	120 kg

Allow +/- 20mm tolerances