

Staebel Haven Rise Lift Chair

Product Code: 5572

The ultra-comfortable, multi-position electric lift chair facilitating access to care and transfers.



The lift function of the Staebel Haven Rise Lift Chair is very appreciated by people who have difficulty getting up, but who also find this function a real help when sitting down. Thanks to its coatings that meet fire-retardant standards and are easy to clean, the HavenRise Chair can be used both at home and in the community. It consists of removable and washable zipped quilt modules at 30°.

The Haven Rise Lift Chair is equipped with removable and / or retractable armrests to facilitate seated – standing transfers (evolving chair in the event of pathological development).

Product Features:

- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°.
- Soft and ultra-comfortable duvets filled with siliconized hollow fiber
- 4 positions: sitting – relax (TV) – rest – lifter
- Different biocompatible coverings available: PVC/PU with M2 fire retardant treatment – mottled fabric – water-repellent velvet
- Different colours variable available depending on the coating
- Allows you to stand and walk tall, instead of leaning forward
- Safe with superior manoeuvrability
- Adjustable footrests and handle heights
- French manufacturing and design.
- One or Two Motors Available
- guarantees: chair: 2 years – cylinder and mechanics: 5 years
- Accessories Available including: mobility kit – patient lift access kit – side shelf – magazine holder – High Density foam – head protection – battery kit

Specifications:

Product Code	5572-1M	5572-2M	5572-1M-XXL
Model	One Motor	Two Motor	XXL
Seat Width	440 / 550 mm	440 / 550 mm	650 mm
Seat Depth	450 / 500 mm	450 / 500 mm	550 mm
Seat Height	520 mm	520 mm	520 mm
Overall Width	930 mm	930 mm	1040 mm
Overall Height	1100 mm	1100 mm	1100 mm
Chair Weight	59 kg	61 kg	76 kg
SWL	160 kg	160 kg	220 kg

