

Staebel Haven Lux Lift Chair

Product Code: 5573

The Staebel Haven Lux Lift Chair combines a compact size with exceptional comfort for home or aged care environments.



The Staebel Haven Lux Lift Chair combines a compact size with exceptional comfort, featuring a low floor-to-seat height of 47 cm for effortless sitting and standing. This elegant, ultra-comfortable chair offers multiple positioning options and has a soft, silky feel due to its 100% polyester velvet covering. With removable, machine-washable zipped duvet modules, it is easy to maintain and can be repaired at a reduced cost.

Product Features:

- One or Two Motors available
- 100% polyester velvet covering gives you a soft and silky feel
- HR 35 (High resilience) foam seat provides excellent comfort and well-being
- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°
- Ultra comfortable and soft duvets filled with siliconised hollow fiber
- 4 positions: sitting – relaxing (TV) – resting – lifting
- Various colours available: Amber, Ruby, Sapphire, Galena – 100% polyester / PVC Pebble Gray, Lagoon Blue & Almond Beige
- 2 adjustable seat depths: 46 or 53 cm
- Warranty: chair: 2 years – cylinder and mechanics: 5 years

Specifications:

Product Code	5573-1M	5573-2M
Model	One Motor	Two Motor
Seat Width	440 mm	440 mm
Seat Depth	460 – 530 mm	460 – 530 mm
Seat Height	470 mm	470 mm
Overall Width	800 mm	800 mm
Overall Height	1120 mm	1120 mm
Chair Weight	54 kg	56 kg
SWL	130 kg	130 kg

