

Staebel Serenity Electric Recliner

Product Code: 5574

Designed for people with reduced mobility assisted by a carer with the use of the electric lift function.



The Staebel Serenity Electric Recline is for indoor use only. It has been designed for people with reduced mobility (total or partial) assisted by a carer. The Serenity Electric provides assistance with standing transfer thanks to its electric LIFT function. Its standing transfer assistance (lift function) allows you to sit down or stand up without effort and it allows the user to maintain their autonomy thanks to its remote control. The Serenity Electric has a leg rest whose deployment is linked to the inclination of the chair. This leg rest is aligned with the seat for great comfort. This push chair allows easy access to dining tables, thanks to the height of its armrests.

Product Features:

- 30° electric tilt: synchronisation of seat, backrest, legrest
- Standing transfer aid by raising the seat
- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°
- Soft & ultra comfortable duvets filled with siliconised hollow fiber
- Bio compatible coatings available:
PVC/PU treated with M2 fire retardant – mottled fabric – water-repellent velvet, microfiber velvet
- Various colours available
(variable depending on the coating)
- Warranty: chair: 2 years – cylinder: 5 years



Specifications:

Product Code	5574-380	5574-440	5574-500	5574-560	5574-650
Seat Width	380 mm	440 mm	500 mm	560 mm	650 mm
Seat Depth	480 mm	480 mm	480 mm	480 mm	480 mm
Overall Width	670 mm	670 mm	730 mm	790 mm	880 mm
Seat Height	480 mm	480 mm	480 mm	480 mm	480 mm
Overall Height	1310 mm	1310 mm	1310 mm	1310 mm	1310 mm
Chair Weight	42 kg	42 kg	42 kg	42 kg	42 kg
SWL	130 kg	130 kg	130 kg	130 kg	130 kg