

Staebel Haven Relax Lift Chair

Product Code: 5577

The ultra-comfortable, multi-position electric lift chair facilitating access to care and transfers.



The Staebel Haven Relax Lift Chair is an ultra-comfortable electric lift chair, removable and/or retractable multi-position to facilitate seated – standing transfers (evolving chair in the event of pathological development). The lift function of the Haven Relax Lift Chair is very appreciated by people who find it difficult to get up but who also find this function a real help when sitting down. Thanks to its PVC coatings to fire-retardant standards, the HavenGlide can be used both at home and in the community.

Using the lifter is very simple. With a simple press of the remote control, the chair returns to a seated position and the chair begins its gentle ascent. You can release the pressure on the remote control at any time, your chair stops. The chair is completely safe!

Product Features:

- It consists of removable and washable zipped quilt modules at 30°
- 4 possible positions of seated, relaxed (TV), rest, lifter for standing transfer
- Soft and ultra-comfortable duvets filled with siliconized hollow fiber
- Different biocompatible coverings available: PVC/PU with M2 fire retardant treatment - mottled fabric - water-repellent velvet
- 8 possible colors (variable depending on the coating)
- Different seat widths available
- Options available include: mobility kit - patient lift access kit - side shelf - magazine holder - High Density foam - head protection - battery kit
- Guarantees: chair: 2 years - cylinder and mechanics: 5 years
- French manufacturing and design.

Specifications:

Product Code	5577-1M	5577-2M	5577-1M-XXL
Model	One Motor	Two Motor	XXL
Seat Width	440 / 550 mm	440 / 550 mm	650 mm
Seat Depth	450 / 500 mm	450 / 500 mm	550 mm
Seat Height	520 mm	520 mm	520 mm
Overall Width	730 mm	730 mm	830 mm
Overall Height	1100 mm	1100 mm	1100 mm
Chair Weight	59 kg	61 kg	76 kg
SWL	160 kg	160 kg	220 kg

