

Staebl Serena Manual Recliner

Product Code: 5582

Designed for versatile use, both indoors and outdoors, catering specifically to those with reduced mobility with the assistance from a caregiver.



The Staebel Serena Manual Recliner is designed for versatile use, both indoors and outdoors, catering specifically to people with reduced mobility who are assisted by a caregiver. Its wheels enable smooth movement across different terrains. In the semi-reclined position, the leg rest aligns with the seat, providing optimal comfort; this is achieved through a synchronized mechanism that adjusts the leg rest as the chair reclines. The chair includes a wide push bar, allowing the caregiver to easily manoeuvre it, and is equipped with drum brakes accessible directly from the push bar for convenient control. The Serena Recliner allows easy access to dining tables, thanks to the armrest height.

Product Features:

- Manual synchronisation of seat, backrest and legrest
- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°
- Ultra comfortable duvets filled with siliconized hollow fiber
- Drum brake system on both rear wheels with manual control
- 3 biocompatible coatings available: PVC/PU treated with M2 fire retardant - mottled fabric - water-repellent velvet
- Various colours available (variable depending on the coating)
- Warranty: Chair: 2 years



Specifications:

Product Code	5582-380	5582-440	5582-500	5582-560	5582-650
Seat Width	380 mm	440 mm	500 mm	560 mm	650 mm
Seat Depth	480 mm	480 mm	480 mm	480 mm	480 mm
Overall Width	670 mm	670 mm	730 mm	790 mm	880 mm
Seat Height	520 mm	520 mm	520 mm	520 mm	520 mm
Overall Height	1350 mm	1350 mm	1350 mm	1350 mm	1350 mm
Chair Weight	37.5 kg	37.5 kg	37.5 kg	37.5 kg	37.5 kg
SWL	130 kg	130 kg	130 kg	130 kg	130 kg