

Staebel Serena Electric Recliner

Product Code: 5585

Designed for versatile use, both indoors and outdoors, catering specifically to those with reduced mobility with the assistance from a caregiver.



The Staebel Serena Electric Recliner is designed for versatile use, both indoors and outdoors, catering specifically to people with reduced mobility who are assisted by a caregiver. Its wheels enable smooth movement across different terrains. In the semi-reclined position, the leg rest aligns with the seat, providing optimal comfort; this is achieved through a synchronised mechanism that adjusts the leg rest as the chair reclines. The chair includes a wide push bar, allowing the caregiver to easily manoeuvre it, and is equipped with drum brakes accessible directly from the push bar for convenient control. The Serena Recliner allows easy access to dining tables, thanks to the armrest height.

Product Features:

- 30° electric tilt: synchronisation of seat, backrest, legrest
- Manual synchronisation of seat, backrest and legrest
- Optimal hygiene: replaceable, removable and machine washable duvet modules at 30°
- Ultra comfortable duvets filled with siliconized hollow fiber
- Drum brake system: both rear wheels with manual control
- 3 biocompatible coatings available: PVC/PU treated with M2 fire retardant, mottled fabric, water-repellent velvet
- Various colours available (variable depending on the coating)
- Warranty: Chair: 2 years, cylinder: 5 years



Specifications:

| Product Code | 5583-380 | 5583-440 | 5583-500 | 5583-560 | 5583-650 |
|----------------|----------|----------|----------|----------|----------|
| Seat Width | 380 mm | 440 mm | 500 mm | 560 mm | 650 mm |
| Seat Depth | 480 mm | 480 mm | 480 mm | 480 mm | 480 mm |
| Overall Width | 670 mm | 670 mm | 730 mm | 790 mm | 880 mm |
| Seat Height | 520 mm | 520 mm | 520 mm | 520 mm | 520 mm |
| Overall Height | 1350 mm | 1350 mm | 1350 mm | 1350 mm | 1350 mm |
| Chair Weight | 38 kg | 38 kg | 38 kg | 38 kg | 38 kg |
| SWL | 130 kg | 130 kg | 130 kg | 130 kg | 130 kg |