

Staebel Serena Electric Recliner

Product Code: 5583

Designed for versatile use, both indoors and outdoors, catering specifically to those with reduced mobility with the assistance from a caregiver.



Specifications:

Product Code	5583-380	5583-440	5583-500	5583-560	5583-650
Seat Width	380 mm	440 mm	500 mm	560 mm	650 mm
Seat Depth	480 mm	480 mm	480 mm	480 mm	480 mm
Overall Width	670 mm	670 mm	730 mm	790 mm	880 mm
Seat Height	520 mm	520 mm	520 mm	520 mm	520 mm
Overall Height	1350 mm	1350 mm	1350 mm	1350 mm	1350 mm
Chair Weight	38 kg	38 kg	38 kg	38 kg	38 kg
SWL	130 kg	130 kg	130 kg	130 kg	130 kg

Issue Date: 150425
Revision Date: 280526 | V2

The Staebel Serena Electric Recliner is designed for versatile use, both indoors and outdoors, catering specifically to people with reduced mobility who are assisted by a caregiver. Its wheels enable smooth movement across different terrains. In the semi-reclined position, the leg rest aligns with the seat, providing optimal comfort; this is achieved through a synchronised mechanism that adjusts the leg rest as the chair reclines.

The chair includes a wide push bar, allowing the caregiver to easily manoeuvre it, and is equipped with drum brakes accessible directly from the push bar for convenient control. The Serena Recliner allows easy access to dining tables, thanks to the armrest height.

Product Features:

- It features a leg rest that extends when the chair is tilted. In the semi-reclined position, the leg rest is aligned with the seat for greater comfort.
- The width of the leg rest and that of the seat are identical in order to guarantee better support
- Extremely maneuverable, it features a wide push bar that allows the caregiver to maneuver the chair with ease.
- The brake is accessible to the assistant directly at the level of the push bar (drum brakes).
- The remote control allows the user to manage the tilt of their chair independently.
- The footrests can be folded down for greater comfort in an inclined position.
- 30° electric tilt function: synchronisation of seat, backrest, legrest
- Manual synchronisation of seat, backrest and legrest
- Optimal hygiene: replaceable, removable and machine washable duvet modules
- Ultra comfortable duvets filled with siliconized hollow fiber
- Drum brake system: both rear wheels with manual control
- 3 biocompatible coatings available: PVC/PU treated with M2 fire retardant, mottled fabric, water-repellent velvet
- Various colours available (variable depending on the coating)
- Warranty: Chair: 2 years, cylinder: 5 years

Accessories

- **Battery Backup (Code: 5598):**
Allows the chair to be used without being plugged into mains power.

