

# Staebel Embark Utilite ExtraCare Day Chair

Product Code: 5911

The Staebel ExtraCare Utilite is a heavy duty, multi purpose adjustable seat height and angle chair, ideal for use in an area where adjustable seating is required.





The Staebel Embark Extra Care Utilite is a heavy duty, multi purpose adjustable seat height and angle chair. Ideal for use in an area where adjustable seating is required, for example, lounge, dining room, bedroom, activity room, etc. This low back chair provides comfort, stable and practicality at an affordable price.

#### Product Benefits:

- Shallow curved low backrest
- Padded armrests
- Comfortable, durable upholstered seat, back and armrest
- Powdercoated steel frame
- 200kg SWL
- Multiple vinyl colours are available – contact the team for more information

#### Specifications:

Product Code	5911
Seat Width	550 mm
Seat Depth	430 mm
Seat Height	450 – 630 mm
Backrest Height	450 mm
Arm Height	230 mm
Overall Width	670 mm
Overall Depth	420 mm
Overall Height	850 – 1030 mm
SWL	200 kg