

Staebel Embark Superior Day Chair

Product Code: 5933

The Staebel Embark Superior Day Chair features a robust 200 kg capacity, low-resilient foam with 2-way stretch upholstery, and convenient 6 cm tilt-and-move wheels for lasting durability.



The Staebel Embark Superior Day Chair was designed and developed by a team of clinicians. It has a strong frame and a high weight capacity of 200 kg. The chair has low resilient foam and 2-way stretch upholstery on seating surfaces. For increased longevity in demanding environments it has 6cm tilt and move wheels and other small improvements.

Product Benefits:

- Low resilient foam and 2-way stretch upholstery
- Padded armrests
- Adjustable seat height and depth
- Powder-coated steel frame
- High backrest incorporating lumbar and lateral support
- Low back option also available.

Specifications:

Product Code	5933
Seat Width	500 mm
Seat Depth	450 - 530 mm
Seat Height	450 - 600 mm
Backrest Height	790 mm
Arm Height	230 mm
Overall Width	700 mm
Overall Depth	730 - 800 mm
Overall Height	1200 - 1350 mm
SWL	200 kg

Issue Date: 281124
Revision Date: 110626 | V3