

Kitchen Perch with Arms

Product Code: 5941

Designed with an angled padded seat which allows the user to rest and then stand again with minimum effort.





An attractive perching stool with an angled padded seat which allow the user to rest and then stand again with minimum effort - particularly helpful with those who tire easily. Suitable for use in the kitchen, bathroom and any activity where a perched seating position is of assistance. The height of the seat can be adjusted by altering the length of the legs.

Product Benefits:

- Shallow curved low backrest
- Armrests as part of the frame
- Comfortable vinyl upholstered seat and back
- Adjustable seat height
- Powdercoated steel frame

Specifications:

Product Code	5941
Seat Width	390 mm
Seat Depth	320 mm
Seat Height	540 - 690 mm
Backrest Height	390 mm
Overall Width	570 mm
SWL	100 kg