

# Staebel Empire High Back Day Chair

Product Code: 5979

The Staebel Empire High Back Day Chair is designed for comfort, mobility, style, and safety for the user, while providing practical solutions for carers.





With adjustable height and seat depth, along with Tilt 'N' Move wheels, it ensures ease of use and flexibility. Its user- and carer-friendly operation adds to its appeal, and as an Australian-made and owned product, it reflects quality and durable craftsmanship.

### Product Features:

- Adjustable Height and Seat depth
- Tilt 'N' Move Wheels
- Black Powdercoated Steel Frame for durability
- User and Carer friendly operation
- Customisable Frame and Padded Section
- Safe Working Load 350kg
- WARRANTY: 1 year - Upholstery, 2 year - Frame

### Specifications:

<b>Product Code</b>	5979-700	5979-800
<b>Width Between Arms</b>	690 mm	790 mm
<b>Seat Depth</b>	450 - 510 mm	450 - 510 mm
<b>Seat Height</b>	440 - 560 mm	440 - 560 mm
<b>Armrest Height</b>	190 mm	190 mm
<b>Backrest Height</b>	790 mm	790 mm
<b>Overall Width</b>	930 mm	1030 mm
<b>Overall Height</b>	1230 mm	1230 mm
<b>Overall Depth</b>	820 mm	820 mm
<b>SWL</b>	350 kg	350 kg