

Staebel Empire High Back Day Chair

Product Code: 5979

The Staebel Empire High Back Day Chair is designed for comfort, mobility, style, and safety for the user, while providing practical solutions for carers.





With adjustable height and seat depth, along with Tilt 'N' Move wheels, it ensures ease of use and flexibility. Its user- and carer-friendly operation adds to its appeal, and as an Australian-made and owned product, it reflects quality and durable craftsmanship.

Product Features:

- Adjustable Height and Seat depth
- Tilt 'N' Move Wheels
- Black Powdercoated Steel Frame for durability
- User and Carer friendly operation
- Customisable Frame and Padded Section
- Safe Working Load 350kg
- WARRANTY: 1 year - Upholstery, 2 year - Frame

Specifications:

Product Code	5979-700	5979-800
Width Between Arms	690 mm	790 mm
Seat Depth	450 - 510 mm	450 - 510 mm
Seat Height	440 - 560 mm	440 - 560 mm
Armrest Height	190 mm	190 mm
Backrest Height	790 mm	790 mm
Overall Width	930 mm	1030 mm
Overall Height	1230 mm	1230 mm
Overall Depth	820 mm	820 mm
SWL	350 kg	350 kg

Issue Date: 290525
Revision Date: 100925 | V2