

Staebel Embark Deluxe Wide Day Chair

Product Code: 5931

The Staebel Embark Deluxe Wide Day Chair has an impressive 250kg weight capacity and extra wide 60cm width, making it ideal for the heavier client.



The Staebel Embark Deluxe Wide Day Chair has an impressive 250kg weight capacity and extra wide 60cm width, making it ideal for the heavier client. It still provides all the benefits of the Embark Deluxe, e.g. adjustable seat depth and height, and comes standard with tilt and move wheels.

Product Benefits:

- High backrest incorporating lumbar and lateral support
- Padded armrests
- Comfortable vinyl upholstered seat and back
- Adjustable seat height and depth
- Powdercoated steel frame with reinforcement for increased load
- A range of vinyl colours available

Specifications:

Product Code	5931
Seat Width	600 mm
Seat Depth	450 - 530 mm
Seat Height	470 - 620 mm
Backrest Height	765 mm
Arm Height	220 mm
Overall Width	800 mm
Overall Depth	740 - 820 mm
Overall Height	1180 - 1330 mm
SWL	250 kg

Issue Date: 281124
Revision Date: 110626 | V3