

Staebel Embark Superior Day Chair

Product Code: 5935

The Staebel Embark Superior Day Chair was designed and developed by a team of clinicians. It has a strong frame and a high weight capacity of 200kg.



The Staebel Embark Superior Day Chair was designed and developed by a team of clinicians. It has a strong frame and a high weight capacity of 200 kg. The chair has low resilient foam and vinyl upholstery on seating surfaces. For increased longevity in demanding environments it has 6cm tilt and move wheels and other small improvements.

Product Benefits:

- Low back option also available.
- Comfortable vinyl seat and back
- Padded armrests
- Adjustable seat height and depth
- Powder-coated steel frame
- High backrest incorporating lumbar and lateral support
- Available in a range of vinyl colours.

Specifications:

Product Code	5935
Seat Width	500 mm
Seat Depth	450 - 530 mm
Seat Height	450 - 600 mm
Backrest Height	790 mm
Arm Height	230 mm
Overall Width	700 mm
Overall Depth	730 - 800 mm
Overall Height	1200 - 1350 mm
SWL	200 kg

Issue Date: 040425
Revision Date: 110626 | V2