

Horizon 3004 Change Table

Product Code: 7-3004

The Horizon 3004 Change Table is ideal for children up to 90kg. The foot control operation provides flexible and easy to control height adjustment, while the single tower lifting column conceals electrics within the column for safe and smooth functioning.



Product Features:

- Electric Hi-Low
- Foot control operation
- 125mm braking castors
- Charcoal Padded Top
- APO Grey & Ironstone Powdercoat Frame

Accessories: *(Added after manufacture)*

- Vertical Side Rails
- Padded Cotside Protectors
- Lap Belt

Specifications:

Product Code	7-3004
Overall Width	700 mm
Overall Length	1400 mm
Hi-Lo Range	480 - 880 mm
SWL	90 kg

A change table for smaller users or children in all healthcare environments.

Optional Features – Added at time of manufacture

- Plas Coating Upgrade
 - Grey (Enhances durability & corrosion resistance of the metal surface)
- Hand Control
- Battery Backup
- Drop Down Padded Side Rails
- Custom Vinyl Colour Padded Top
- Custom Colour Powdercoated Frame

Powdercoat Colours Available

Ironstone	APO Grey	Rivergum Beige	Satin Black	Precious Silver Pearl
Pearl White	Signal Red	Orange X15	Lemon Yellow	Pink Blush
Dark Violet	French Blue	Navy	Cabana Green	Mint Green

Vinyl Colours Available

Black	Amethyst	Smoke	Pacific
Teal	Cloud	Frost	Kiwi
Poppy	Orange	Lemon	Snow

