

Horizon 3250 Change Table

Product Code: 7-3250

The Horizon 3250 Change Table is ideal for multiple environments. Designed & manufactured in Australia, our 3250 Change Table has a vertical lift function and a variety of optional features available including a backrest, central locking castors & foot control.



Suitable for
Healthcare Environments



SWL
250kg



Australian
Made & Owned



Customisable
Designs

The Horizon 3250 Change Table is a versatile and heavy-duty care solution designed to support safe, comfortable, and efficient changing, treatment, and examination procedures across a wide range of healthcare and community environments. With a safe working load of 250kg, it is ideal for facilities requiring dependable performance for children, adults, and bariatric care applications.

Designed and manufactured in Australia, the 3250 features smooth vertical lift height adjustment to assist carers in achieving an ergonomic working position, helping reduce manual handling strain while improving accessibility for users. The stable lifting platform and generously padded top provide comfort, security, and confidence during daily care routines.

Built for demanding environments, the Horizon 3250 combines a durable powder-coated steel frame with high-quality castors for easy manoeuvrability and secure positioning when in use. Its practical design makes it well suited to hospitals, disability services, schools, therapy spaces, aged care facilities, and specialist hygiene rooms.

Product Features:

- Australian Made & Owned
- Charcoal Vinyl Padded Top
- APO Grey & Ironstone powdercoat frame
- 125mm braking and directional castors
- Hand Control
- Safe Working Load of 250kg

Optional Features - Added at time of manufacture

- Plas Coating Upgrade - Grey (Enhances durability & corrosion resistance of the metal surface)
- Backrest
- Battery Backup
- Jumbo Battery System
- Central Locking Castors
- Drop Down Padded Side Rails
- Face Hole
- Foot Control
- Custom Colour Vinyl Top
- Custom Colour Powdercoated Frame

Accessories: (Added after manufacture)

- Vertical Rails
- Cotside Protectors
- Lap Belt



Vertical Rails
(7-380VRCT-G)



Lap Belt
(7-LAP/B)



Drop Down Padded Side Rails



Custom Vinyl Padded Top & Powdercoat Frame

Specifications:

Style	Without Backrest	With Backrest
Height Range	450 - 900 mm	450 - 900 mm
Overall Width	700 mm	700 mm
Overall Length	1880 mm	1880 mm
Backrest Angle	-	75°
SWL	250 kg	250 kg

Issue Date: 280125
Revision Date: 080626 | V3