

Horizon 4000 Patient Trolley

Produce Code: 7-4PT

Designed for use as both a day surgery bed and a recovery trolley. The Trendelenburg permits the carer to regulate the blood circulation of the patient while being transferred from operating room to ward.



Product Features:

- Australian Made & Owned
- Electric Hi Low
- Electric Backrest with Emergency CPR Release
- Electric Forward & Reverse Trendelenburg
- 150mm Twin wheel anti-static Central Locking Castors
- Push Handles
- Battery Backup
- Corner Buffers
- APO Grey & Ironstone powdercoat frame
- Warranty:
10 years frame, 5 years electrics, 1 year castors &
1 year battery

Specifications:

Overall Dimensions	840 x 2230 mm
Sleeping Surface	750 x 2025 mm
Height Range	390 - 830 mm
Underbed Clearance	160 mm
Back Angle	75°
TR / ATR	16° / 16°
Trolley Weight	120 kg
SWL	300 kg



Optional Features – Added at time of manufacture

- Kneebreak
- Attendant Control
- Jumbo battery system
- Underbed Light (Warm White)
- Underbed Light with Custom Colours
- USB Charger Ports
- 200mm castors upgrade (increases min/max height)

Accessories – Added after manufacture

- 400HR Horizontal Rails
- Oxygen Bottle Carrier
- Self Help Pole & Triangle
- Vertical Rails
- I.V. Pole