

Horizon Montar Wall Mounted Change Table

Product Code: 7-MONTAR

The Montar Wall Mounted Change Table offers space saving efficiencies and will integrate seamlessly into your facility. It includes state of the art electrics to achieve a 385mm low height and an 985mm high height.



Product Features:


- Australian Made & Owned
- Charcoal Vinyl Padded Top
- Corded Handset
- Drop Down Safety Rail Folds Up For Saving Space
- Electric Hi-Low
- Wall Mounted
- APO Grey Powdercoating
- Battery Back Up
- IP rating of IP X6
- Warranty: 10 years frame, 5 years electrics, 1 year castors & 1 year battery

Accessories: *(Added after manufacture)*

- Lap belt

Specifications:

Size	Paediatric	Adult
Sleeping Surface	700 x 1400 mm	700 x 1880 mm
Overall Width (Open / Closed)	350 / 1030 mm	350 / 1030 mm
Overall Length	1400 mm	1880 mm
Rail Height (From Padded Top)	110 mm	110 mm
Height Range (Approximate)	385 - 985 mm	385 - 985 mm
Product Weight	100 kg	100 kg
SWL	200 kg	200 kg



A space saving wall mounted change table for all care environments.

Optional Features - Added at time of manufacture

- Custom Width & Length
- Bluetooth Handset
- Custom Vinyl Colour
- Custom Powdercoat Colours

Vinyl Colours Available:

Black	Amethyst	Smoke	Pacific
Teal	Cloud	Frost	Kiwi
Poppy	Orange	Lemon	Snow

