

# Horizon Montar Splash Wall Mounted Change Table

Product Code: 7-MONTSS

The Montar Splash Wall Mounted Change Table offers space saving efficiencies and will integrate seamlessly into your facility. It includes state of the art electrics and is compatible for wet areas with a durable, stainless steel construction.



## Product Features:

- Australian Made & Owned
- Navy Vinyl Padded Top
- Corded Handset
- Navy Powder coated Drop Down Safety Rails that fold up for saving space
- Electric Hi-Low
- Wall Mounted
- APO Powder Coated Stainless Steel Frame
- Ti Motion Electrics
- Battery Back Up
- IP Rating of IP66

## Accessories: *(Added after manufacture)*

- Lap belt

## Specifications:

Size	Paediatric	Adult
<b>Sleeping Surface</b>	700 x 1400 mm	700 x 1880 mm
<b>Overall Width (Open / Closed)</b>	350 / 1030 mm	350 / 1030 mm
<b>Overall Length</b>	1400 mm	1880 mm
<b>Rail Height (From Padded Top)</b>	110 mm	110 mm
<b>Height Range (Approximate)</b>	385 - 985 mm	385 - 985 mm
<b>Product Weight</b>	100 kg	100 kg
<b>SWL</b>	200 kg	200 kg



Water  
Repellent



SWL  
200kg



Stainless Steel  
Construction



Customisable  
Designs

# A space saving wall mounted change table for all care environments.

## Optional Features - Added at time of manufacture

- Custom Width & Length
- Bluetooth Handset
- Custom Vinyl Colour
- Custom Powdercoat Colours

## Vinyl Colours Available:

Black	Amethyst	Smoke	Pacific
Teal	Cloud	Frost	Kiwi
Poppy	Orange	Lemon	Snow

