

Horizon Montar Splash Wall Mounted Change Table

Product Code: 7-MONTSS

The Montar Splash Wall Mounted Change Table offers space saving efficiencies and will integrate seamlessly into your facility. It includes state of the art electrics and is compatible for wet areas with a durable, stainless steel construction.



*All images shown are of the Adult Model
(1880mm Length)



Product Features:

- Australian Made & Owned
- Navy Vinyl Padded Top
- Corded Handset
- Navy Powder coated Drop Down Safety Rails that fold up for saving space
- Electric Hi-Low
- Wall Mounted
- APO Powder Coated Stainless Steel Frame
- Ti Motion Electrics
- Battery Back Up
- IP Rating of IP66

Accessories: (Added after manufacture)

- Lap belt



Custom Vinyl Colours



Safety Rail

Issue Date: 140425
Revision Date: 190226 | V2

Specifications:

Size	Paediatric	Adult
Sleeping Surface	700 x 1400 mm	700 x 1880 mm
Overall Width (Open / Closed)	350 / 1030 mm	350 / 1030 mm
Overall Length	1400 mm	1880 mm
Rail Height (From Padded Top)	110 mm	110 mm
Height Range (Approximate)	385 - 985 mm	385 - 985 mm
Bed Base	Padded Vinyl	Padded Vinyl
Under Bed Base	APO Grey Powdercoated Steel & Vinyl	Plastic Cover
Product Weight	100 kg	100 kg
SWL	200 kg	200 kg

Optional Features - Added at time of manufacture

- Custom Width & Length
- Custom Powdercoat Colours
- Custom Vinyl Colour
- Bluetooth Handset
- Safety Rail



Water
Repellent



SWL
200kg



Stainless Steel
Construction



Customisable
Designs