

# Staebel Gusto Wheelchair

Product Code: 46GBW

Gusto Wheelchair caters for the bariatric client designed for comfort, and is ideal in the aged care and hospital environment.



The Staebel Gusto Heavy Duty Wheelchair caters for the bariatric client. Featuring padded vinyl upholstery, heavy duty folding frame, 24 inch rear wheels and inset foot plates for client comfort, this wheelchair is ideal in the aged care and hospital environment. It comes standard with urethane tyres, eliminating the need to maintain tyre pressure. The swing away leg rests are removable and also adjust in length. The armrests will remove to enable side transfers if necessary.

## Product Features:

- Footplates are alloy for longevity and ease of cleaning.
- Removable, swingaway, adjustable length legrests
- Removable armrests
- Padded vinyl upholstery
- 24" rear mag type wheels, urethane tyres

## Specifications:

<b>Product Code</b>	46108	46608	47108
<b>Seat Width</b>	<b>610 mm</b>	<b>660 mm</b>	<b>710 mm</b>
<b>Seat Depth</b>	460 mm	520 mm	550 mm
<b>Seat Height</b>	520 mm	520 mm	520 mm
<b>Backrest Height</b>	430 mm	480 mm	520 mm
<b>Armrest Height</b>	220 mm	235 mm	220 mm
<b>Overall Width</b>	850 mm	880 mm	930 mm
<b>Overall Depth</b>	1100 mm	1150 mm	1170 mm
<b>Overall Height</b>	930 mm	1000 mm	990 mm
<b>Fold Width</b>	350 mm	350 mm	350 mm
<b>Net/Gross Weight</b>	26 / 30 kg	36 / 43 kg	36 / 46 kg
<b>SWL</b>	240 kg	270 kg	290 kg



240 - 290kg  
SWL



Bariatric  
Needs