

Secura 2 Mental Health Bed

Product Code: 7-SEC2

Providing dignity and comfort for mental health patients, with a safe working load of 300kg, designed specifically to reduce stress on carers by relieving physical strain. Suited for low to medium care patients who may have clinical care needs with reduced anti-ligature features.



*Image for illustration purposes. The castors on the Secura 2 come with Twin wheel 125mm castor.



Product Features:

- Electric Hi-Low
- Enclosed anti-ligature area
- All bolts secured with Loctite
- Global Lock Out with Fob
- Embedded footboard control
- Concealed wiring
- Braided electrical power supply
- Secure bolt on head and footboards
- 125mm Twin Wheel Braking & Directional Castors
- Australian owned and Made



Enclosed
Mattress Supports



Anti Tamper
Screws



Natural Oak
Curved Endboards



Ironstone
Powdercoat

Specifications:

Size	Single	King Single
Overall Dimensions	920 x 2070 mm	1090 x 2070 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm
Hi-Lo Position	355 - 805 mm	355 - 805 mm
Underbed Clearance	135 mm	135 mm
SWL	300 kg	300 kg



Backrest is optional

Ensuring high comfort & safety in mental health environments.

Optional Features – Added at time of manufacture


- Backrest
- Knee break
- Trendelenberg
- Battery Back up
- Key power lockout
- CPR Release
- 200mm Fixed Build in Extension
- Custom powdercoat colours
- Custom Endboard colours
- Central Locking (with mechanism covers)
- Removable Brake Lever
- Removable Mains Cable

Endboard Colours & Styles


Custom Endboard, styles and colours also available.



Laminate Curved Style



Natural Oak




Oyster Grey

Accessories – Added after manufacture

- Mattresses
- Mobile Battery Supply
- Bluetooth Handset (7-BT-4-HB200)



Central Locking
(with mechanism covers)



Key Power Lockout



CPR Release



Removable Brake Lever



Custom Designs

Issue Date: 050224
Revision Date: 260226 | V3