

Staebel Horizontal 360HR Bed Rails

Product Code: 7-360HR

Designed to minimise the risk of falls, which are still (unfortunately) a leading cause of injuries in hospitals, nursing homes, and private homes where aging people are being cared for.



Product Features:

- They require no special tools to install.
- Sold as a pair
- Ironstone powder coat
- Secured with a spring loaded lock nut
- 155mm between railings
- Suitable for Staebel P5500 Homecare Bed

Specifications:

Product Code	7-360HR
Height	400 mm
Length	1400 mm